

Antibodies Product Datasheet

Product Name:	Cofilin Mouse mAb(Mix-mA)
Catalog Number:	N1498
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cofilin
Concentration:	1 mg/ml
Gene Name:	CFL1 CFL
Protein Name:	Cofilin
Other Name:	Cofilin-1 (18 kDa phosphoprotein) (p18) (Cofilin, non-muscle isoform)
Human Swiss Prot No.:	P23528
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.71% sodium azide.
Specificity:	This antibody detects endogenous levels of Cofilin at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year