
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | GAS41 Antibody |
| Catalog Number: | C10237 |
| Reactivity: | Human, Mouse |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human GAS41. AA range:1-50 |
| Concentration: | 1 mg/ml |
| Gene Name: | YEATS4 |
| Protein Name: | YEATS domain-containing protein 4 |
| Other Name: | YEATS4; GAS41; YEATS domain-containing protein 4; Glioma-amplified sequence 41; Gas41; NuMA-binding protein 1; NuBI-1; NuBI1 |
| Human Swiss Prot No.: | O95619 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | GAS41 Polyclonal Antibody detects endogenous levels of GAS41 protein. |
| Storage and Stability: | -20°C/1 year |