

Antibodies Product Datasheet

Product Name:	ERO1L Antibody
Catalog Number:	C48687
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ERO1L
Concentration:	1 mg/ml
Gene Name:	ERO1L UNQ434/PRO865
Protein Name:	ERO1-like protein alpha (ERO1-L) (ERO1-L-alpha) (EC 1.8.4.-) (Endoplasmic oxidoreductin-1-like protein) (Oxidoreductin-1-L-alpha)
Other Name:	
Human Swiss Prot No.:	Q96HE7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ERO1L at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year