

Antibodies Product Datasheet

Product Name:	RFA2 (Phospho-Ser33) Antibody
Catalog Number:	P12-1172
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	RFA2 (Phospho-Ser33) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human RFA2 around the phosphorylation site of Ser33.
Concentration:	1 mg/ml
Gene Name:	RFA2
Other Name:	REPA2; replication factor-A protein 2; replication protein A 32 kDa subunit; replication protein A2; RF-A; RP-A; RPA2; RPA32
Human Swiss Prot No.:	P15927
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	RFA2 (Phospho-Ser33) Antibody detects endogenous levels of RFA2 only when phosphorylated at Ser33.
Storage and Stability:	Stable at -20°C for at least 1 year.