

## Antibodies Product Datasheet

<b>Product Name:</b>	RIP2 Antibody
<b>Catalog Number:</b>	B8180
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human RIPK2. AA range:146-195
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	RIPK2
<b>Protein Name:</b>	Receptor-interacting serine/threonine-protein kinase 2
<b>Other Name:</b>	RIPK2; CARDIAK; RICK; RIP2; Receptor-interacting serine/threonine-protein kinase 2; CARD-containing interleukin-1 beta-converting enzyme-associated kinase; CARD-containing IL-1 beta ICE-kinase; RIP-like-interacting CLARP kinase; Receptor-in
<b>Human Swiss Prot No.:</b>	O43353
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	RIP2 Polyclonal Antibody detects endogenous levels of RIP2 protein.
<b>Storage and Stability:</b>	-20°C/1 year