

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | NRBP2 Antibody |
| Catalog Number: | C45073 |
| Reactivity: | Human, Mouse |
| Applications: | WB, ELISA |
| Dilutions: | WB 1:500-2000 ELISA 1:5000-20000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human protein . at AA range: 40-120 |
| Concentration: | 1 mg/ml |
| Gene Name: | NRBP2 PP9320 TRG16 |
| Protein Name: | Nuclear receptor-binding protein 2 (Transformation-related gene 16 protein) (TRG-16) |
| Other Name: | |
| Human Swiss Prot No.: | Q9NSY0 |
| Formulation: | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Specificity: | NRBP2 Polyclonal Antibody detects endogenous levels of protein. |
| Storage and Stability: | -20°C/1 year |