

Antibodies Product Datasheet

Product Name:	V1RL2 Antibody
Catalog Number:	G792
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VN1R2. AA range:88-137
Concentration:	1 mg/ml
Gene Name:	VN1R2
Protein Name:	Vomeroneasal type-1 receptor 2
Other Name:	VN1R2; V1RL2; Vomeroneasal type-1 receptor 2; G-protein coupled receptor GPCR25; hGPCR25; V1r-like receptor 2
Human Swiss Prot No.:	Q8NFZ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	V1RL2 Polyclonal Antibody detects endogenous levels of V1RL2 protein.
Storage and Stability:	-20°C/1 year