

Antibodies Product Datasheet

Product Name:	Anti-Human hCG β (ABT-HCG β) Mouse Monoclonal Antibody
Catalog Number:	V0081
Clone ID:	ABT-HCG β
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CGB3/CGB5/CGB8/CGB7/CGB1/CGB2
Gene Name:	P0DN86/P0DN87/A6NKQ9/Q6NT52
Other Name:	CG beta, CG-beta, CgB, CGB3, CGB3_HUMAN, CGB5, CGB7, CGB7_HUMAN, CGB8, CGHB_HUMAN, Choriogonadotropin subunit beta 3, Choriogonadotropin subunit beta 7, Choriogonadotropin subunit beta, Chorionic gonadotropin beta polypeptide, Chorionic gonadotropin chain beta
Isotype	IgG1, Kappa
Localization:	Placenta
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize β subunits of human CG, and has no cross reaction with LH, FSH or TSH proteins.
Storage and Stability:	-20°C/1 year