

Antibodies Product Datasheet

Product Name:	PMR1 Monoclonal Antibody
Catalog Number:	N0523
Reactivity:	Human,Monkey
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PMR1 expressed in E. Coli.
Concentration:	
Gene Name:	ATP2C1
Protein Name:	Calcium-transporting ATPase type 2C member 1
Other Name:	ATP2C1; KIAA1347; PMR1L; HUSSY-28; Calcium-transporting ATPase type 2C member 1; ATPase 2C1; ATP-dependent Ca(2+) pump PMR1
Human Swiss Prot No.:	P98194
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PMR1 Monoclonal Antibody detects endogenous levels of PMR1 protein.
Storage and Stability:	-20°C/1 year