

Antibodies Product Datasheet

Product Name:	VATL Antibody
Catalog Number:	C51163
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal VATL
Concentration:	1 mg/ml
Gene Name:	ATP6V0C ATP6C ATP6L ATPL
Protein Name:	V-type proton ATPase 16 kDa proteolipid subunit (V-ATPase 16 kDa proteolipid subunit) (Vacuolar proton pump 16 kDa proteolipid subunit)
Other Name:	V-type proton ATPase 16 kDa proteolipid subunit (V-ATPase 16 kDa proteolipid subunit) (Vacuolar proton pump 16 kDa proteolipid subunit)
Human Swiss Prot No.:	P27449
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of VATL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year