

Antibodies Product Datasheet

Product Name:	Cadherin-20 Antibody
Catalog Number:	C12095
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDH20. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	CDH20
Protein Name:	Cadherin-20
Other Name:	CDH20; CDH7L3; Cadherin-20
Human Swiss Prot No.:	Q9HBT6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cadherin-20 Polyclonal Antibody detects endogenous levels of Cadherin-20 protein.
Storage and Stability:	-20°C/1 year