

Antibodies Product Datasheet

Product Name:	MAN2A1 Antibody
Catalog Number:	C49815
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MAN2A1
Concentration:	1 mg/ml
Gene Name:	MAN2A1 MANA2
Protein Name:	Alpha-mannosidase 2 (EC 3.2.1.114) (Golgi alpha-mannosidase II) (AMan II) (Man II) (Mannosidase alpha class 2A member 1) (Mannosyl-oligosaccharide 1,3-1,6-alpha-mannosidase)
Other Name:	
Human Swiss Prot No.:	Q16706
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAN2A1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year