

## Antibodies Product Datasheet

<b>Product Name:</b>	NFkB p65 Monoclonal Antibody(5G6)
<b>Catalog Number:</b>	N1376
<b>Reactivity:</b>	Human,Mouse,Rat
<b>Applications:</b>	WB,IHC-p,IF/ICCIP
<b>Dilutions:</b>	WB: 1:1000-3000 IP:1:200 IF 1:200 IHC 1:50-300
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Transcription factor p65
<b>Concentration:</b>	
<b>Gene Name:</b>	RELA
<b>Protein Name:</b>	Transcription factor p65
<b>Other Name:</b>	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
<b>Human Swiss Prot No.:</b>	Q04206
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous p65 proteins.
<b>Storage and Stability:</b>	-20°C/1 year