

Antibodies Product Datasheet

Product Name:	DGK- δ Antibody
Catalog Number:	C11140
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKD. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	DGKD
Protein Name:	Diacylglycerol kinase delta
Other Name:	DGKD; KIAA0145; Diacylglycerol kinase delta; DAG kinase delta; 130 kDa diacylglycerol kinase; Diglyceride kinase delta; DGK-delta
Human Swiss Prot No.:	Q16760
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- δ Polyclonal Antibody detects endogenous levels of DGK- δ protein.
Storage and Stability:	-20°C/1 year