

Antibodies Product Datasheet

Product Name:	CREB-1 Antibody
Catalog Number:	B7052
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CREB. AA range:96-145
Concentration:	1 mg/ml
Gene Name:	CREB1
Protein Name:	Cyclic AMP-responsive element-binding protein 1
Other Name:	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1
Human Swiss Prot No.:	P16220
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CREB-1 Polyclonal Antibody detects endogenous levels of CREB-1 protein.
Storage and Stability:	-20°C/1 year