
Antibodies Product Datasheet

Product Name:	EYFP-Tag Antibody
Catalog Number:	C40006
Reactivity:	Species independent
Applications:	WB
Dilutions:	WB: 1:1000 -3000
Modification:	Unmodified/Total
Source:	Polyclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of EYFP-Tag
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects EYFP and YFP tag fusion proteins.
Storage and Stability:	-20°C/1 year