

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | DDR1 Monoclonal Antibody |
| Catalog Number: | N0193 |
| Reactivity: | Human |
| Applications: | WB,ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant fragment of DDR1 (aa602-681) expressed in E. Coli. |
| Concentration: | |
| Gene Name: | DDR1 |
| Protein Name: | Epithelial discoidin domain-containing receptor 1 |
| Other Name: | DDR1; CAK; EDDR1; NEP; NTRK4; PTK3A; RTK6; TRKE; Epithelial discoidin domain-containing receptor 1; Epithelial discoidin domain receptor 1; CD167 antigen-like family member A; Cell adhesion kinase; Discoidin receptor tyrosine kinase; HGK2; |
| Human Swiss Prot No.: | Q08345 |
| Formulation: | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol. |
| Specificity: | DDR1 Monoclonal Antibody detects endogenous levels of DDR1 protein. |
| Storage and Stability: | -20°C/1 year |