

Antibodies Product Datasheet

Product Name:	β -actin Monoclonal Antibody(5B7), HRP Conjugated
Catalog Number:	N1151
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Insect
Applications:	IF/ICC,WB,IHC-p
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:2000-5000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG1
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	ACTB
Protein Name:	Actin cytoplasmic 1
Other Name:	ACTB
Human Swiss Prot No.:	P60709
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	β -actin Monoclonal Antibody(5B7) specially designed for your Western blot as internal control.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.