

Antibodies Product Datasheet

Product Name:	ERK 3 Antibody
Catalog Number:	C11136
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p97 MAPK. AA range:461-510
Concentration:	1 mg/ml
Gene Name:	MAPK6
Protein Name:	Mitogen-activated protein kinase 6
Other Name:	MAPK6; ERK3; PRKM6; Mitogen-activated protein kinase 6; MAP kinase 6; MAPK 6; Extracellular signal-regulated kinase 3; ERK-3; MAP kinase isoform p97; p97-MAPK
Human Swiss Prot No.:	Q16659
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERK 3 Polyclonal Antibody detects endogenous levels of ERK 3 protein.
Storage and Stability:	-20°C/1 year