

Antibodies Product Datasheet

Product Name:	UBR1 Antibody
Catalog Number:	R12-3665
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	UBR1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human UBR1.
Concentration:	1 mg/ml
Gene Name:	UBR1
Other Name:	E3 ubiquitin-protein ligase UBR1; N-recogin-1; Ubiquitin-protein ligase E3-alpha-1; Ubiquitin-protein ligase E3-alpha-I; UBR1
Human Swiss Prot No.:	Q8IWW7
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	UBR1 Antibody antibody detects endogenous levels of UBR1.
Storage and Stability:	Stable at -20°C for at least 1 year.