

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | MORF4L2 Antibody |
| Catalog Number: | C50378 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit,IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human MORF4L2 |
| Concentration: | 1 mg/ml |
| Gene Name: | MORF4L2 KIAA0026 MRGX |
| Protein Name: | Mortality factor 4-like protein 2 (MORF-related gene X protein) (Protein MSL3-2) (Transcription factor-like protein MRGX) |
| Other Name: | |
| Human Swiss Prot No.: | Q15014 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of MORF4L2 at Human, Mouse,Rat |
| Storage and Stability: | Store at -20° C/1 year |