

Antibodies Product Datasheet

Product Name:	DPP9 Antibody
Catalog Number:	C51063
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DPP9
Concentration:	1 mg/ml
Gene Name:	DPP9 DPRP2
Protein Name:	Dipeptidyl peptidase 9 (DP9) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 2) (DPRP-2) (Dipeptidyl peptidase IX) (DPP IX) (Dipeptidyl peptidase-like protein 9) (DPLP9)
Other Name:	
Human Swiss Prot No.:	Q86TI2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DPP9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year