
Antibodies Product Datasheet

Product Name:	IPMK Antibody
Catalog Number:	C11376
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IPMK. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	IPMK
Protein Name:	Inositol polyphosphate multikinase
Other Name:	IPMK; IMPK; Inositol polyphosphate multikinase; Inositol 1; 3,4,6-tetrakisphosphate 5-kinase
Human Swiss Prot No.:	Q8NFU5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IPMK Polyclonal Antibody detects endogenous levels of IPMK protein.
Storage and Stability:	-20°C/1 year