

Antibodies Product Datasheet

Product Name:	gp210 Antibody
Catalog Number:	C50158
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human gp210
Concentration:	1 mg/ml
Gene Name:	NUP210 KIAA0906 PSEC0245
Protein Name:	Nuclear pore membrane glycoprotein 210 (Nuclear pore protein gp210) (Nuclear envelope pore membrane protein POM 210) (POM210) (Nucleoporin Nup210) (Pore membrane protein of 210 kDa)
Other Name:	
Human Swiss Prot No.:	Q8TEM1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of gp210 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year