

Antibodies Product Datasheet

Product Name:	FGFR3 (Phospho-Tyr760) Antibody
Catalog Number:	A8103
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGFR3 (Phospho-Tyr760)
Concentration:	1 mg/ml
Gene Name:	FGFR3 JTK4
Protein Name:	FGFR3 (Phospho-Tyr760)
Other Name:	Fibroblast growth factor receptor 3 (FGFR-3) (EC 2.7.10.1) (CD antigen CD333)
Human Swiss Prot No.:	P22607
Formulation:	Liquid in PBS containing 50% glycerol, and 0.104% sodium azide.
Specificity:	This antibody detects endogenous levels of FGFR3 (Phospho-Tyr760) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year