
Antibodies Product Datasheet

Product Name:	c-Abl Monoclonal Antibody
Catalog Number:	N0086
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of c-Abl (aa577-650) expressed in E. Coli.
Concentration:	
Gene Name:	ABL1
Protein Name:	Tyrosine-protein kinase ABL1
Other Name:	ABL1; ABL; JTK7; Tyrosine-protein kinase ABL1; Abelson murine leukemia viral oncogene homolog 1; Abelson tyrosine-protein kinase 1; Proto-oncogene c-Abl; p150
Human Swiss Prot No.:	P00519
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	c-Abl Monoclonal Antibody detects endogenous levels of c-Abl protein.
Storage and Stability:	-20°C/1 year