

Antibodies Product Datasheet

Product Name:	HP1- β (Phospho-Thr51) Antibody
Catalog Number:	A9091
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HP1- β (Phospho-Thr51)
Concentration:	1 mg/ml
Gene Name:	CBX1 CBX
Protein Name:	HP1- β (Phospho-Thr51)
Other Name:	Chromobox protein homolog 1 (HP1Hsbeta) (Heterochromatin protein 1 homolog beta) (HP1 beta) (Heterochromatin protein p25) (M31) (Modifier 1 protein) (p25beta)
Human Swiss Prot No.:	P83916
Formulation:	Liquid in PBS containing 50% glycerol, and 0.185% sodium azide.
Specificity:	This antibody detects endogenous levels of HP1- β (Phospho-Thr51) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year