

## Antibodies Product Datasheet

<b>Product Name:</b>	CD22 Antibody
<b>Catalog Number:</b>	B8253
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000, ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human CD22. at AA range: 780-860
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CD22 SIGLEC2
<b>Protein Name:</b>	CD22
<b>Other Name:</b>	B-cell receptor CD22 (B-lymphocyte cell adhesion molecule) (BL-CAM) (Sialic acid-binding Ig-like lectin 2) (Siglec-2) (T-cell surface antigen Leu-14) (CD antigen CD22)
<b>Human Swiss Prot No.:</b>	P20273
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of CD22
<b>Storage and Stability:</b>	-20°C/1 year