

Antibodies Product Datasheet

Product Name:	Gemin3 Antibody
Catalog Number:	C17863
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DDX20. AA range:273-322
Concentration:	1 mg/ml
Gene Name:	DDX20
Protein Name:	Probable ATP-dependent RNA helicase DDX20
Other Name:	DDX20; DP103; GEMIN3; Probable ATP-dependent RNA helicase DDX20; Component of gems 3; DEAD box protein 20; DEAD box protein DP 103; Gemin-3
Human Swiss Prot No.:	Q9UHI6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gemin3 Polyclonal Antibody detects endogenous levels of Gemin3 protein.
Storage and Stability:	-20°C/1 year