

## Antibodies Product Datasheet

<b>Product Name:</b>	Adducin $\alpha/\beta$ Antibody
<b>Catalog Number:</b>	B0002
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IP, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human ADD1. AA range:688-737
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ADD1/ADD2
<b>Protein Name:</b>	Alpha-adducin/Beta-adducin
<b>Other Name:</b>	ADD1; ADDA; Alpha-adducin; Erythrocyte adducin subunit alpha; ADD2; ADDB; Beta-adducin; Erythrocyte adducin subunit beta
<b>Human Swiss Prot No.:</b>	P35611/P35612
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Adducin $\alpha/\beta$ Polyclonal Antibody detects endogenous levels of Adducin $\alpha/\beta$ protein.
<b>Storage and Stability:</b>	-20°C/1 year