

Antibodies Product Datasheet

Product Name:	TFIIIB90-1 Antibody
Catalog Number:	C11452
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TF3B. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	BRF1
Protein Name:	Transcription factor IIIB 90 kDa subunit
Other Name:	BRF1; BRF; GTF3B; TAF3B2; TAF3C; Transcription factor IIIB 90 kDa subunit; TFIIIB90; hTFIIIB90; B-related factor 1; BRF-1; hBRF; TAF3B2; TATA box-binding protein-associated factor; RNA polymerase III, subunit 2
Human Swiss Prot No.:	Q92994
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TFIIIB90-1 Polyclonal Antibody detects endogenous levels of TFIIIB90-1 protein.
Storage and Stability:	-20°C/1 year