

Antibodies Product Datasheet

Product Name:	Anti-Human Melan-A (ABT193) Mouse Monoclonal Antibody
Catalog Number:	V0393
Clone ID:	ABT193
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	MLANA
Gene Name:	Q16655
Other Name:	Antigen LB39 AA, Antigen LB39-AA, Antigen SK29 AA, Antigen SK29-AA, MAR1_HUMAN, MART 1, MART-1, MART1, Melan A, Melan A protein, Melanoma antigen recognized by T cells 1, Melanoma antigen recognized by T-cells 1, MLAN A, MLANA, OTTHUMP00000021036, OTTHUMP00000021037, OTTHUMP00000021038, Protein Melan-A
Isotype	IgG1, Kappa
Localization:	Skin/ Melanoma
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Melan-A protein.
Storage and Stability:	-20°C/1 year