
Antibodies Product Datasheet

Product Name:	VEGF-A Antibody
Catalog Number:	R12-2407
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VEGF-A. AA range:110-159
Concentration:	1 mg/ml
Gene Name:	VEGFA
Protein Name:	Vascular endothelial growth factor A
Other Name:	VEGFA; VEGF; Vascular endothelial growth factor A; VEGF-A; Vascular permeability factor; VPF
Human Swiss Prot No.:	P15692
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VEGF-A Polyclonal Antibody detects endogenous levels of VEGF-A protein.
Storage and Stability:	-20°C/1 year