
Antibodies Product Datasheet

Product Name:	SCAMP1 Antibody
Catalog Number:	C18492
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SCAMP1. AA range:273-322
Concentration:	1 mg/ml
Gene Name:	SCAMP1
Protein Name:	Secretory carrier-associated membrane protein 1
Other Name:	SCAMP1; SCAMP; Secretory carrier-associated membrane protein 1; Secretory carrier membrane protein 1
Human Swiss Prot No.:	O15126
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SCAMP1 Polyclonal Antibody detects endogenous levels of SCAMP1 protein.
Storage and Stability:	-20°C/1 year