

Antibodies Product Datasheet

Product Name:	Anti-Human MCM2 (ABT-MCM2) Mouse Monoclonal Antibody
Catalog Number:	V0093
Clone ID:	ABT-MCM2
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	MCM2
Gene Name:	P49736
Other Name:	BM28, CCNL 1, CCNL1, CDC like 1, cdc19, CDCL 1, CDCL1, Cell division cycle like 1, Cyclin like 1, D3S3194, DNA replication licensing factor MCM2, KIAA0030, MCM 2, MCM2, MCM2 minichromosome maintenance deficient 2 mitotin, MCM2 minichromosome maintenance deficient 2 mitotin (S. cerevisiae), MCM2 minichromosome maintenance deficient 2, mitotin, MCM2_HUMAN, MGC10606, Minichromosome maintenance complex component 2, Minichromosome maintenance deficient 2 (mitotin), Minichromosome maintenance deficient 2 mitotin, Minichromosome maintenance protein 2, Minichromosome maintenance protein 2 homolog, Mitotin, Nuclear protein BM28, OTTHUMP00000216047, OTTHUMP00000216050
Isotype	IgG1, Kappa
Localization:	Tonsil
Positive Tissue Control:	Nuclear
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human MCM2 protein. In western blotting of HEK293 cell lysate, the antibody can label a 130 kDa band.
Storage and Stability:	-20°C/1 year