
Antibodies Product Datasheet

Product Name:	Rb (Phospho-Ser807/811) Antibody
Catalog Number:	A9795
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Rb (Ser807 and 811)
Concentration:	1 mg/ml
Gene Name:	RB1
Protein Name:	Rb (Ser807/811)
Other Name:	Retinoblastoma-associated protein (p105-Rb) (pRb) (Rb) (pp110)
Human Swiss Prot No.:	P06400
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Rb (phospho-Ser807 or 811)
Storage and Stability:	-20°C/1 year