
Antibodies Product Datasheet

Product Name:	GLP-1R Antibody
Catalog Number:	G273
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GLP1R. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	GLP1R
Protein Name:	Glucagon-like peptide 1 receptor
Other Name:	GLP1R; Glucagon-like peptide 1 receptor; GLP-1 receptor; GLP-1-R; GLP-1R
Human Swiss Prot No.:	P43220
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GLP-1R Polyclonal Antibody detects endogenous levels of GLP-1R protein.
Storage and Stability:	-20°C/1 year