

Antibodies Product Datasheet

Product Name:	Wee 1 (Phospho-Ser53) Antibody
Catalog Number:	A8199
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WEE1 around the phosphorylation site of Ser53. AA range:19-68
Concentration:	1 mg/ml
Gene Name:	WEE1
Protein Name:	Wee1-like protein kinase
Other Name:	WEE1; Wee1-like protein kinase; WEE1hu; Wee1A kinase
Human Swiss Prot No.:	P30291
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Wee 1 (S53) Polyclonal Antibody detects endogenous levels of Wee 1 protein only when phosphorylated at S53.
Storage and Stability:	-20°C/1 year