

Antibodies Product Datasheet

Product Name:	XRCC1 (Phospho-Thr284) Antibody
Catalog Number:	A9130
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human XRCC1 (Phospho-Thr284)
Concentration:	1 mg/ml
Gene Name:	XRCC1
Protein Name:	XRCC1 (Phospho-Thr284)
Other Name:	DNA repair protein XRCC1 (X-ray repair cross-complementing protein 1)
Human Swiss Prot No.:	P18887
Formulation:	Liquid in PBS containing 50% glycerol, and 0.171% sodium azide.
Specificity:	This antibody detects endogenous levels of XRCC1 (Phospho-Thr284) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year