

Antibodies Product Datasheet

Product Name:	FGR (Phospho-Tyr412) Antibody
Catalog Number:	A9807
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGR (Phospho-Tyr412)
Concentration:	1 mg/ml
Gene Name:	FGR SRC2
Protein Name:	FGR (Phospho-Tyr412)
Other Name:	Tyrosine-protein kinase Fgr (EC 2.7.10.2) (Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog) (Proto-oncogene c-Fgr) (p55-Fgr) (p58-Fgr) (p58c-Fgr)
Human Swiss Prot No.:	P09769
Formulation:	Liquid in PBS containing 50% glycerol, and 0.107% sodium azide.
Specificity:	This antibody detects endogenous levels of FGR (Phospho-Tyr412) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year