

## Antibodies Product Datasheet

<b>Product Name:</b>	TIP60 Antibody
<b>Catalog Number:</b>	B0814
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human TIP60. AA range:52-101
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	KAT5
<b>Protein Name:</b>	Histone acetyltransferase KAT5
<b>Other Name:</b>	KAT5; HTATIP; TIP60; Histone acetyltransferase KAT5; 60 kDa Tat-interactive protein; Tip60; Histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; Lysine acetyltransferase 5; cPLA(2)-interacting protein
<b>Human Swiss Prot No.:</b>	Q92993
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	TIP60 Polyclonal Antibody detects endogenous levels of TIP60 protein.
<b>Storage and Stability:</b>	-20°C/1 year