

Antibodies Product Datasheet

Product Name:	CRF-RII Antibody
Catalog Number:	G073
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CRHR2. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	CRHR2
Protein Name:	Corticotropin-releasing factor receptor 2
Other Name:	CRHR2; CRF2R; CRH2R; Corticotropin-releasing factor receptor 2; CRF-R-2; CRF-R2; CRFR-2; Corticotropin-releasing hormone receptor 2; CRH-R-2; CRH-R2
Human Swiss Prot No.:	Q13324
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRF-RII Polyclonal Antibody detects endogenous levels of CRF-RII protein.
Storage and Stability:	-20°C/1 year