

Antibodies Product Datasheet

Product Name:	MRI1 Antibody
Catalog Number:	C50920
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MRI1
Concentration:	1 mg/ml
Gene Name:	MRI1 MRDI UNQ6390/PRO21135
Protein Name:	Methylthioribose-1-phosphate isomerase (M1Pi) (MTR-1-P isomerase) (EC 5.3.1.23) (Mediator of RhoA-dependent invasion) (S-methyl-5-thioribose-1-phosphate isomerase) (Translation initiation factor eIF-2B subunit alpha/beta/delta-like protein)
Other Name:	
Human Swiss Prot No.:	Q9BV20
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MRI1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year