

Antibodies Product Datasheet

Product Name:	G6PC2 Antibody
Catalog Number:	C43854
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human G6PC2
Concentration:	1 mg/ml
Gene Name:	G6PC2 IGRP
Protein Name:	G6PC2
Other Name:	Glucose-6-phosphatase 2 (G-6-Pase 2) (G6Pase 2) (EC 3.1.3.9) (Islet-specific glucose-6-phosphatase catalytic subunit-related protein)
Human Swiss Prot No.:	Q9NQR9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.67% sodium azide.
Specificity:	This antibody detects endogenous levels of G6PC2 at Human/Mouse
Storage and Stability:	-20°C/1 year