

Antibodies Product Datasheet

Product Name:	GATA-1 (Acetyl-Lys312) Antibody
Catalog Number:	D12264
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GATA-1 (Acetyl Lys312)
Concentration:	1 mg/ml
Gene Name:	GATA1 ERYF1 GF1
Protein Name:	GATA-1 (Acetyl Lys312)
Other Name:	Erythroid transcription factor (Eryf1) (GATA-binding factor 1) (GATA-1) (GF-1) (NF-E1 DNA-binding protein)
Human Swiss Prot No.:	P15976
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat GATA-1 (Acetyl Lys312)
Storage and Stability:	-20°C/1 year