

## Antibodies Product Datasheet

<b>Product Name:</b>	CARD8 Antibody
<b>Catalog Number:</b>	C49762
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human CARD8
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CARD8 KIAA0955 NDPP1
<b>Protein Name:</b>	Caspase recruitment domain-containing protein 8 (Apoptotic protein NDPP1) (CARD-inhibitor of NF-kappa-B-activating ligand) (CARDINAL) (DACAR) (Tumor up-regulated CARD-containing antagonist of CASP9) (TUCAN)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9Y2G2
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of CARD8 at Human
<b>Storage and Stability:</b>	Store at -20° C/1 year