

Antibodies Product Datasheet

Product Name:	MSL3 Antibody
Catalog Number:	C50619
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MSL3
Concentration:	1 mg/ml
Gene Name:	MSL3 MSL3L1
Protein Name:	Male-specific lethal 3 homolog (Male-specific lethal-3 homolog 1) (Male-specific lethal-3 protein-like 1) (MSL3-like 1)
Other Name:	
Human Swiss Prot No.:	Q8N5Y2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MSL3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year