

Antibodies Product Datasheet

Product Name:	POLR2B Antibody
Catalog Number:	C50844
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human POLR2B
Concentration:	1 mg/ml
Gene Name:	POLR2B
Protein Name:	DNA-directed RNA polymerase II subunit RPB2 (EC 2.7.7.6) (DNA-directed RNA polymerase II 140 kDa polypeptide) (DNA-directed RNA polymerase II subunit B) (RNA polymerase II subunit 2) (RNA polymerase II subunit B2)
Other Name:	
Human Swiss Prot No.:	P30876
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POLR2B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year