

Antibodies Product Datasheet

Product Name:	JMJD6(N-term) mouse mAb
Catalog Number:	N0769
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human JMJD6(N-terminus) fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	jmjd6
Protein Name:	
Other Name:	Apoptotic cell clearance receptor; Bifunctional arginine demethylase and lysyl-hydroxylase JMJD6; Histone arginine demethylase JMJD6; JmjC domain-containing protein 6; JMJD 6; JMJD6; JMJD6_HUMAN; Jumonji domain containing 6; Jumonji domain-containing protein 6; Kl
Human Swiss Prot No.:	Q6NYC1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of JMJD6(N-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year