

## Antibodies Product Datasheet

<b>Product Name:</b>	CD292 Antibody
<b>Catalog Number:</b>	C30487
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from the N-terminal region of human BMPR1A. AA range:1-50
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	BMPR1A
<b>Protein Name:</b>	Bone morphogenetic protein receptor type-1A
<b>Other Name:</b>	BMPR1A; ACVRLK3; ALK3; Bone morphogenetic protein receptor type-1A; BMP type-1A receptor; BMPR-1A; Activin receptor-like kinase 3; ALK-3; Serine/threonine-protein kinase receptor R5; SKR5; CD292
<b>Human Swiss Prot No.:</b>	P36894
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	CD292 Polyclonal Antibody detects endogenous levels of CD292 protein.
<b>Storage and Stability:</b>	-20°C/1 year